



SHARPER IMAGE®

AIR COMPRESSION BOOTS

Please read all instructions and warnings prior to use. Keep this manual as it contains important information for future reference.


WHAT YOU GET

- a. Air Compression Boots
- b. Control Unit
- c. Air Hose Extension
- d. Power Adapter with DC Jack


a.




b.



c.



d.



• Accessory type and quantity may vary depending on the model. •

• The images shown in this manual are for illustrative purposes only and may be subject to change. •


THE FUTURE OF RECOVERY

01

Plug the barrel jack on the end of the power cord into the input port on the control unit.

Plug the power adapter into a wall outlet.

Note: Do not use the power adapter to charge or power other devices.



02

Unzip the boot or boots you want to use for your massage. The boots are labeled L (left leg) and R (right leg).

Find a comfortable place to sit or recline with your legs extended for the duration of your massage and remove any footwear (wearing socks is recommended). Next, pull each boot over the corresponding leg and then zip the boot up all the way to the top of the zipper track.

Lock the zippers in place by flipping the zipper tabs down fully flat against the zipper track. The tabs should point towards the toes of the boots. Failure to properly lock the zippers may cause the boot to open during inflation, reducing compression.

IMPORTANT: Never attempt to stand or walk while wearing the boots. The boots are intended for massage only. They are not designed for stability or mobility. Using the boots for purposes other than their intended use may damage the boots or result in personal injury.

03

IMPORTANT: Make sure the hose release button is facing upwards when connecting each hose. Press each connector into each connection port slowly and firmly until it clicks securely into place.

To use both boots at the same time, plug the end of the air hose extension with a single connector into the air hose connection port on the control unit. Next, connect both boot air hoses to the connection ports on the other end of the air hose extension.


To use a single boot, connect the boot air hose directly to the control unit or connect it using the air hose extension. When using the air hose extension with a single boot, seal the unused connection port with the rubber plug located on the bottom of the connector junction (the top side features the Sharper Image logo).

The air hoses can be disconnected by pressing each hose release button while gently pulling the connector away from the connection port. Always grip the plastic connectors directly. Never pull or push using the hoses.


Air Hose Extension (Top View)



Air Hose Extension (Bottom View)



Single-Boot Usage



Both-Boots Usage

To power on the massager, press and hold the power button. To power off the massager, press and hold the power button again.

Quick press the timer button after compression starts, pressing the timer button to increase the massage timers (add 10mins, 20mins, return to initial time).

- Press the mode button to cycle through the 5 massage modes.
- Pulse (default) – all zones compress and then decompress
 - Gradual – zones each compress sequentially and then all zones decompress
 - Sequence – zones each compress and decompress sequentially
 - Zones – select up to 5 zones for targeted pulse massage
 - Manual – select up to 5 zones and manually adjust compression.


In zones mode or manual mode, make your compression area choices using the up, down and select buttons.

After all mode selections have been made, press and hold the select button to start your massage.

Quick press the compression button to cycle through the 5 compression levels.

Quick press the power button during a massage to begin decompression. The massager will automatically power off after 3 minutes and the boots will be fully decompressed.

To stop the decompression, press and hold the power button to turn off the device.



NOTE: The images shown here are for illustrative purposes only and may be subject to change, the actual descriptions of control panel and their locations may vary depending on the model.

CARE INSTRUCTIONS

- Spot clean only. Unplug and wipe with a damp cloth.
- Air-dry only. Never place in a dryer.
- Allow to dry fully before connecting to power sources.
- Do not submerge in water.
- Do not use detergent or soap when cleaning.
- Hang the boots or lay them flat for long-term storage.
- Never crumple or wad the boots.
- Do not crimp or forcefully bend any of the included hoses.

SAFETY WARNINGS

PLEASE READ BEFORE USE

- Not a children's item. Keep out of the reach of children.
- As with any massage, bruising may occur. Test your comfort level and check targeted areas often. Stop your massage immediately if you feel any increase in physical discomfort or pain.
- If you experience numbness, tingling, or loss of sensation, stop your massage immediately. Seek professional medical advice if these symptoms persist.
- Do not fall asleep while using the massager.
- Do not use attachments of any kind other than those that may be included.
- Do not leave this massager unattended when it is plugged into a power source.
- Do not expose to extreme temperatures or moisture.
- Keep away from sources of water and heated surfaces.
- Never touch the power cable or ports with wet hands.
- Do not carry this item by the power cable, hoses, or attachments.
- For indoor use only.
- Thoroughly inspect product before every use. If the product appears worn, frayed, splintered, cracked, or broken in any manner, discontinue use and discard immediately.
- This product is not intended to treat, cure, or prevent medical conditions.
- Do not use this item if you are pregnant or if have an injury or illness. Consult your physician if you have a medical condition that may prevent safe use of this item.

WARNING: Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

FCC STATEMENT

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Warning: Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

SHARPER IMAGE® NAME AND LOGOS (INCLUDING THE 77 LOGO, TOMORROW'S TOMORROW, AND POWERBOOST) ARE TRADEMARKS OR REGISTERED TRADEMARKS. DISTRIBUTED BY AND © 2022 MERCHSOURCE, LLC. ALL RIGHTS RESERVED.

PRINTED IN MEXICO
SC:12-16-2024 / Elbrd: J.M
OBPN:1018195
SHNMDL: 1018195_B
ADD: 1018196

8 1 4 1 0 0 3 0 6

